

ERBB4 Antibody



PACO09124

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ERBB4

Uniprot

Q15303

Synonyms:

v-erb-a erythroblastic leukemia viral oncogene homolog 4 (avian); ERBB4; HER4; MGC138404; p180erbB4 ;

Immunogen:

Human ERBB4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A