

ERCC6L Antibody



PACO09132

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ERCC6L

Uniprot

Q2NWX8

Synonyms:

excision repair cross-complementing rodent repair deficiency, complementation group 6-like; ERCC6L; FLJ20105; MGC131695; PICH ;

Immunogen:

Human ERCC6L

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A