

# ERO1L Antibody



PACO09140

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ERO1A

**Uniprot**

Q96HE7

**Synonyms:**

ERO1-like (*S. cerevisiae*); ERO1L; ERO1-alpha ;

**Immunogen:**

Human ERO1L.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A