

# ERO1LB Antibody



PACO09141

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ERO1B

**Uniprot**

Q86YB8

**Synonyms:**

ERO1-like beta (S. cerevisiae);

**Immunogen:**

Human ERO1LB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A