

ESRRG Antibody



PACO09148

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ESRRG

Uniprot

P62508

Synonyms:

estrogen-related receptor gamma; ESRRG; DKFZp781L1617; ERR3; FLJ16023; KIAA0832; NR3B3 ; ERRgamma

Immunogen:

Human ESRRG.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A