

ETFA Antibody



PACO09151

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ETFA

Uniprot

P13804

Synonyms:

electron-transfer-flavoprotein, alpha polypeptide; ETFA; EMA; GA2; MADD ;

Immunogen:

Human ETFA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A