

EVI2A Antibody



PACO09159

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

EVI2A

Uniprot

P22794

Synonyms:

ecotropic viral integration site 2A; EVI2A; EVDA; EVI2 ;

Immunogen:

Human EVI2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A