

EXOC1 Antibody



PACO09162

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

EXOC1

Uniprot

Q9NV70

Synonyms:

exocyst complex component 1; EXOC1; BM-102; FLJ10893; SEC3; SEC3L1; SEC3P ;

Immunogen:

Human EXOC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A