

EXOSC10 Antibody



PACO09168

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

EXOSC10

Uniprot

Q01780

Synonyms:

exosome component 10; EXOSC10; PM-Scl; PM/Scl-100; PMSCL; PMSCL2; RRP6; Rrp6p; p2; p3; p4 ;

Immunogen:

Human EXOSC10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A