

EXOSC2 Antibody



PACO09169

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

EXOSC2

Uniprot

Q13868

Synonyms:

exosome component 2; EXOSC2; RRP4; Rrp4p; hRrp4p; p7 ;

Immunogen:

Human EXOSC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A