

EXTL2 Antibody



PACO09175

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

EXTL2

Uniprot

Q9UBQ6

Synonyms:

exostoses (multiple)-like 2; EXTL2; EXTR2 ;

Immunogen:

Human EXTL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A