

F12 Antibody



PACO09181

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

F12

Uniprot

P00748

Synonyms:

coagulation factor XII (Hageman factor); F12; HAE3; HAEX; HAF ;

Immunogen:

Human F12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A