

FABP3 Antibody



PACO09186

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FABP3

Uniprot

P05413

Synonyms:

fatty acid, binding protein 3, muscle and heart (mammary-derived growth inhibitor); FABP3; FABP11; H-FABP; MDGI; O-FABP ;

Immunogen:

Human FABP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A