

FABP4 Antibody



PACO09187

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FABP4

Uniprot

P15090

Synonyms:

fatty acid, binding protein 4, adipocyte; FABP4; A-FABP; aP2 ;

Immunogen:

Human FABP4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A