

# FAM107A Antibody



PACO09198

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAM107A

**Uniprot**

O95990

**Synonyms:**

family with sequence similarity 107, member A; FAM107A; DRR1; FLJ30158; FLJ45473; TU3A ;

**Immunogen:**

Human FAM107A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A