

# FAM111B Antibody



PACO09199

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAM111B

**Uniprot**

Q6SJ93

**Synonyms:**

family with sequence similarity 111, member B; FAM111B; CANP ;

**Immunogen:**

Human FAM111B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A