

METTL21A Antibody



PACO09202

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

METTL21A

Uniprot

Q8WXB1

Synonyms:

family with sequence similarity 119, member A; FAM119A; HCA557b; MGC45373 ;

Immunogen:

Human FAM119A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A