

# FAM129A Antibody



PACO09207

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

FAM129A

**Uniprot**

Q9BZQ8

**Synonyms:**

family with sequence similarity 129, member A; FAM129A; C1orf24; FLJ38228; GIG39; NIBAN ;

**Immunogen:**

Human FAM129A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A