

FAM49B Antibody



PACO09216

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

FAM49B

Uniprot

Q9NUQ9

Synonyms:

family with sequence similarity 49, member B; FAM49B; BM-009; DKFZp686B04128; L1 ;

Immunogen:

Human FAM49B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A