

FAM57B Antibody



PACO09218

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FAM57B

Uniprot

Q71RH2

Synonyms:

family with sequence similarity 57, member B; FAM57B; DKFZp434I2117; FP1188 ;

Immunogen:

Human FAM57B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A