

FAM72A Antibody



PACO09222

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FAM72A

Uniprot

Q5TYM5

Synonyms:

family with sequence similarity 72, member A; FAM72A; MGC46651; MGC57827 ;

Immunogen:

Human FAM72A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A