

FANCD2 Antibody



PACO09228

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FANCD2

Uniprot

Q9BXW9

Synonyms:

Fanconi anemia, complementation group D2; FANCD2; DKFZp762A223; FA-D2; FA4; FACD; FAD; FAD2; FANCD; FLJ23826 ;

Immunogen:

Human FANCD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A