

PACO09239

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FAU

Uniprot

P62861

Synonyms:

Finkel-Biskis-Reilly murine sarcoma virus (FBR-MuSV) ubiquitously expressed; FAU; FAU1; FLJ22986; Fub1; Fubi; MNSFbeta; RPS30; asr1 ;

Immunogen:

Human FAU.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A