

FBXL13 Antibody



PACO09246

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FBXL13

Uniprot

Q8NEE6

Synonyms:

F-box and leucine-rich repeat protein 13; FBXL13; FLJ38068; Fbl13; MGC21636 ;

Immunogen:

Human FBXL13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A