

FBXL18 Antibody



PACO09248

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FBXL18

Uniprot

Q96ME1

Synonyms:

F-box and leucine-rich repeat protein 18; FBXL18; FLJ10776; FLJ11467; FLJ26934; FLJ38075; FLJ41541; Fbl18 ;

Immunogen:

Human FBXL18.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A