

FBXL8 Antibody



PACO09249

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FBXL8

Uniprot

Q96CD0

Synonyms:

F-box and leucine-rich repeat protein 8; FBXL8; FBL8; FLJ11278; MGC19959 ;

Immunogen:

Human FBXL8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A