

PACO09250

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FBXO15

Uniprot

Q8NCQ5

Synonyms:

F-box protein 15; FBXO15; FBX15; MGC39671 ;

Immunogen:

Human FBXO15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A