

FBXO22 Antibody



PACO09253

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FBXO22

Uniprot

Q8NEZ5

Synonyms:

F-box protein 22; FBXO22; FBX22; FISTC1; FLJ13986; MGC31799 ;

Immunogen:

Human FBXO22.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A