

FBXO9 Antibody



PACO09260

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FBXO9

Uniprot

Q9UK97

Synonyms:

F-box protein 9; FBXO9; DKFZp434C0118; FBX9; KIAA0936; NY-REN-57; VCIA1; dJ341E18.2 ;

Immunogen:

Human FBXO9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A