

FGD2 Antibody



PACO09280

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FGD2

Uniprot

Q7Z6J4

Synonyms:

FYVE, RhoGEF and PH domain containing 2; FGD2; FLJ00048; FLJ40929; MGC71330; ZFYVE4 ;

Immunogen:

Human FGD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A