

FGF13 Antibody



PACO09285

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FGF13

Uniprot

Q92913

Synonyms:

fibroblast growth factor 13; FGF13; FGF2; FHF-2; FHF2 ;

Immunogen:

Human FGF13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A