

# FGFBP2 Antibody



PACO09286

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FGFBP2

**Uniprot**

Q9BYJ0

**Synonyms:**

fibroblast growth factor binding protein 2; FGFBP2; HBP17RP; KSP37 ;

**Immunogen:**

Human FGFBP2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A