

# FGFR1OP2 Antibody



PACO09288

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FGFR1OP2

**Uniprot**

Q9NVK5

**Synonyms:**

FGFR1 oncogene partner 2; FGFR1OP2; DKFZp564O1863; HSPC123-like ;

**Immunogen:**

Human FGFR1OP2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A