

FGFR4 Antibody



PACO09290

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FGFR4

Uniprot

P22455

Synonyms:

fibroblast growth factor receptor 4; FGFR4; CD334; JTK2; MGC20292; TKF ;

Immunogen:

Human FGFR4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A