

FHL1 Antibody



PACO09295

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FHL1

Uniprot

Q13642

Synonyms:

four and a half LIM domains 1; Fhl1; KyoT ;

Immunogen:

Human Fhl1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A