

FKBP7 Antibody



PACO09308

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FKBP7

Uniprot

Q9Y680

Synonyms:

FK506 binding protein 7; FKBP7; FKBP23; MGC9420; PPIase ;

Immunogen:

Human FKBP7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A