

FMNL1 Antibody



PACO09320

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FMNL1

Uniprot

O95466

Synonyms:

formin-like 1; FMNL1; C17orf1; C17orf1B; FHOD4; FMNL; KW-13; MGC133052; MGC1894; MGC21878 ;

Immunogen:

Human FMNL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A