

FMR1 Antibody



PACO09326

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

FMR1

Uniprot

Q06787

Synonyms:

fragile X mental retardation 1; FMR1; FMRP; FRAXA; MGC87458; POF; POF1 ;

Immunogen:

Human FMR1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A