

FNTA Antibody



PACO09331

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

FNTA

Uniprot

P49354

Synonyms:

farnesyltransferase, CAAX box, alpha; FNTA; FPTA; MGC99680; PGGT1A; PTAR2 ;

Immunogen:

Human FNTA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A