

FOXE3 Antibody



PACO09335

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

FOXE3

Uniprot

Q13461

Synonyms:

forkhead box E3; FOXE3; ASMD; FKHL12; FREAC8 ;

Immunogen:

Human FOXE3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A