

FRS3 Antibody



PACO09348

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FRS3

Uniprot

O43559

Synonyms:

fibroblast growth factor receptor substrate 3; FRS3; FRS2B; FRS2beta; MGC17167; SNT-2; SNT2 ;

Immunogen:

Human FRS3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A