

# FUS Antibody



PACO09359

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

FUS

**Uniprot**

P35637

**Synonyms:**

fusion (involved in t(12;16) in malignant liposarcoma); FUS; CHOP; FUS-CHOP; FUS1; TLS; TLS/CHOP; hnRNP-P2 ;

**Immunogen:**

Human FUS.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A