

FUT8 Antibody



PACO09363

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FUT8

Uniprot

Q9BYC5

Synonyms:

fucosyltransferase 8 (alpha (1,6) fucosyltransferase); FUT8; MGC26465 ;

Immunogen:

Human FUT8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A