

# GABARAPL2 Antibody



PACO09379

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GABARAPL2

**Uniprot**

P60520

**Synonyms:**

GABA(A) receptor-associated protein-like 2; GABARAPL2; ATG8; GATE-16; GATE16; GEF2 ;

**Immunogen:**

Human GABARAPL2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A