

GABRA3 Antibody



PACO09385

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GABRA3

Uniprot

P34903

Synonyms:

gamma-aminobutyric acid, (GABA) A receptor, alpha 3; GABRA3; MGC33793 ;

Immunogen:

Human GABRA3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A