

# GADD45GIP1 Antibody



PACO09393

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GADD45GIP1

**Uniprot**

Q8TAE8

**Synonyms:**

growth arrest and DNA-damage-inducible, gamma interacting protein 1; GADD45GIP1; CKBBP2; CRIF1; MGC4667; MGC4758; PLINP-1; PRG6; Plinp1 ;

**Immunogen:**

Human GADD45GIP1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A