

GALNT14 Antibody



PACO09402

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GALNT14

Uniprot

Q96FL9

Synonyms:

UDP-N-acetyl-alpha-D-galactosamine: polypeptide N-acetylgalactosaminyltransferase 14 (GalNAc-T14); GALNT14; FLJ12691; FLJ13977; GALNT15; GalNAc-T10; GalNAc-T14 ;

Immunogen:

Human GALNT14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A