

# GBA3 Antibody



PACO09423

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

GBA3

**Uniprot**

Q9H227

**Synonyms:**

glucosidase, beta, acid, 3 (cytosolic); GBA3; CBGL1; GLUC; KLRP; MGC104276; MGC126878 ;

**Immunogen:**

Human GBA3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A