

PACO09431

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GCHFR

**Uniprot**

P30047

**Synonyms:**

GTP cyclohydrolase I feedback regulator; GCHFR; GFRP; HsT16933; MGC138467; MGC138469; P35 ;

**Immunogen:**

Human GCHFR.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A