

PACO09434

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

GCLM

**Uniprot**

P48507

**Synonyms:**

glutamate-cysteine ligase, modifier subunit; GCLM; GLCLR ;

**Immunogen:**

Human GCLM.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A